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1. Abu-Absi, Nicholas R., Brian M. Kenty, Maryann Ehly Cuellar, Michael C. Borys, Sivakesava Sakhamuri, David J. Strachan, Michael C. Hausladen, and Zheng Jian Li. 2011. "Real Time Monitoring of Multiple Parameters in Mammalian Cell Culture Bioreactors Using an In-Line Raman Spectroscopy Probe." *Biotechnology and Bioengineering* 108 (5): 1215–21. <https://doi.org/10.1002/bit.23023>
2. A. Henslee, Erin, Ruth M. Torcal Serrano, Fatima H. Labeed, Rita I. Jabr, Christopher H. Fry, Michael P. Hughes, and Kai F. Hoettges. 2016. "Accurate Quantification of Apoptosis Progression and Toxicity Using a Dielectrophoretic Approach." *Analyst* 141 (23): 6408–15. <https://doi.org/10.1039/C6AN01596D>
3. Ansonge, Sven, Geoffrey Esteban, and Georg Schmid. 2010. "On-Line Monitoring of Responses to Nutrient Feed Additions by Multi-Frequency Permittivity Measurements in Fed-Batch Cultivations of CHO Cells." *Cytotechnology* 62 (2): 121–32. <https://doi.org/10.1007/s10616-010-9267-z>
4. Ansonge, Sven, Olivier Henry, Marc Aucoin, Robert Voyer, John P. Carvell, and Amine Kamen. 2010. "Monitoring the Cell Size Distribution of Mammalian Cell Cultures Using On-Line Capacitance Measurements." In *Cells and Culture*, edited by Thomas Noll, 853–59. Dordrecht: Springer Netherlands. https://doi.org/10.1007/978-90-481-3419-9_150
5. Ansonge, Sven, Stéphane Lanthier, Julia Transfiguración, Olivier Henry, and Amine Kamen. 2011. "Monitoring Lentiviral Vector Production Kinetics Using Online Permittivity Measurements." *Biochemical Engineering Journal* 54 (1): 16–25. <https://doi.org/10.1016/j.bej.2011.01.002>
6. Cadena-Herrera, Daniela, Joshua E. Esparza-De Lara, Nancy D. Ramírez-Ibañez, Carlos A. López-Morales, Néstor O. Pérez, Luis F. Flores-Ortiz, and Emilio Medina-Rivero. 2015. "Validation of Three Viable-Cell Counting Methods: Manual, Semi-Automated, and Automated." *Biotechnology Reports* 7 (September): 9–16. <https://doi.org/10.1016/j.btre.2015.04.004>
7. Dabros, Michal, Danielle Dennewald, David J. Currie, Mark H. Lee, Robert W. Todd, Ian W. Marison, and Urs von Stockar. 2009. "Cole-Cole, Linear and Multivariate Modeling of Capacitance Data for on-Line Monitoring of Biomass." *Bioprocess and Biosystems Engineering* 32 (2): 161–73. <https://doi.org/10.1007/s00449-008-0234-4>
8. Downey, Brandon J., Lisa J. Graham, Jeffrey F. Breit, and Nathaniel K. Glutting. 2014. "A Novel Approach for Using Dielectric Spectroscopy to Predict Viable Cell Volume (VCV) in Early Process Development." *Biotechnology Progress* 30 (2): 479–87. <https://doi.org/10.1002/btpr.1845>
9. Fernandes, Juhi, Jayme Currie, Kevin Ramer, and An Zhang. 2019. "Development of Capacitance Tools: At-Line Method for Assessing Biomass of Mammalian Cell Culture and Fixed Cell Calibration Standard." *Biotechnology Journal* 14 (4): e1800283. <https://doi.org/10.1002/biot.201800283>
10. Justice, C., J. Leber, D. Freimark, P. Pino Grace, M. Kraume, and P. Czermak. 2011. "Online- and Offline- Monitoring of Stem Cell Expansion on Microcarrier." *Cytotechnology* 63 (4): 325–35. <https://doi.org/10.1007/s10616-011-9359-4>
11. Ma, Fuduo, An Zhang, David Chang, Orlin D. Velev, Kelly Wiltberger, and Rashmi Kshirsagar. 2019. "Real-Time Monitoring and Control of CHO Cell Apoptosis by in Situ Multifrequency Scanning Dielectric Spectroscopy." *Process Biochemistry* 80 (May): 138–45. <https://doi.org/10.1016/j.procbio.2019.02.017>
12. Moore, Brandon, Ryan Sanford, and An Zhang. 2019. "Case Study: The Characterization and Implementation of Dielectric Spectroscopy (Biocapacitance) for Process Control in a Commercial GMP CHO Manufacturing Process." *Biotechnology Progress* 35 (3): e2782. <https://doi.org/10.1002/btpr.2782>
13. Pais, Daniel A. M., Chris Brown, Anastasia Neuman, Krishanu Mathur, Ines A. Isidro, Paula M. Alves, and Peter G. Slade. 2020. "Dielectric Spectroscopy to Improve the Production of RAAV Used in Gene Therapy." *Processes* 8 (11): 1456. <https://doi.org/10.3390/pr8111456>
14. Pais, Daniel A. M., Paulo Galrao, Anastasiya Kryzhanska, Jeremie Barbau, Ines A. Isidro, and Paula Alves. 2020. "Holographic Imaging of Insect Cell Cultures: Online Non-Invasive Monitoring of Adeno-Associated Virus Production and Cell Concentration." *Processes* 8 (4): 487. <https://doi.org/10.3390/pr8040487>
15. Pais, Daniel A. M., Rui M. C. Portela, Manuel J. T. Carrondo, Ines A. Isidro, and Paula M. Alves. 2019. "Enabling PAT in Insect Cell Bioprocesses: In Situ Monitoring of Recombinant Adeno-associated Virus Production by Fluorescence Spectroscopy." *Biotechnology and Bioengineering* 116 (11): 2803–14. <https://doi.org/10.1002/bit.27117>
16. Pan, Xiao, Ciska Dalm, Rene H. Wijffels, and Dirk E. Martens. 2017. "Metabolic Characterization of a CHO Cell Size Increase Phase in Fed-Batch Cultures." *Applied Microbiology and Biotechnology* 101 (22): 8101–13. <https://doi.org/10.1007/s00253-017-8531-y>
17. Petiot, Emma, Sven Ansonge, Manuel Rosa-Calatrava, and Amine Kamen. 2017. "Critical Phases of Viral Production Processes Monitored by Capacitance." *Journal of Biotechnology* 242 (January): 19–29. <https://doi.org/10.1016/j.jbiotec.2016.11.010>

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